

## **Talin Full Length Expression and Purification**

pET30a-Talin Full Length  
MW: 271.7KD  
Ext Coef: 0.334  
pI 5.8  
C-terminal His Tag (non cleavable)

### **Expression:**

- 1) 100ml o/n culture (TB/kana).
- 2) Seed 1L TB/kana with o/n culture. Grow at 37°C until OD600 reaches 0.4-0.6. Transfer to 30°C. Induce with 0.2mM IPTG. Grow o/n at 30°C.
- 3) Harvest 8krpm, 10min. Resuspend in His Binding Buffer. LN2 freeze.

### **Purification:**

- 1) Thaw cells (3L). Add 2mM PMSF. Add 1mM b-ME.
- 2) Lyse cells w/lysozyme. Rock 30min, 4°C.
- 3) Add 10mM MgCl<sub>2</sub>. Digest DNA w/DNaseI. Rock 15min, 4°C.
- 4) Add 0.1% Tx-100. Rock 15min, 4°C.
- 5) Spin 18krpm, 1hr, 4°C.
- 6) Wash 3ml Ni-NTA column w/binding buffer.
- 7) Load supernatant onto column.
- 8) Wash column w/50ml binding buffer.
- 9) Wash column w/ 40ml wash buffer (hi salt).
- 10) Elute w/elution buffer.
- 11) Dialize into dialysis buffer.
- 12) Further clean with Hi Load 16/60 Superdex200 Prep Grade column.
- 13) LN2 freeze.

### **Buffers:**

#### **Binding Buffer:**

20mM Tris HCl pH 7.9  
500mM NaCl  
5mM Imidazole  
1mM β-ME

#### **Elution Buffer:**

20mM Tris HCl pH 7.9  
500mM NaCl  
200mM Imidazole  
1mM β-ME

#### **Wash Buffer:**

20mM Tris HCl pH 7.9  
500mM NaCl  
30mM Imidazole  
1mM β-ME

#### **Dialysis Buffer:**

20mM Tris HCl pH 7.4  
150mM NaCl  
1mM EDTA  
5mM β-ME